## **30" Ductless Fume Hood**

## Model # SS-330-DCH

The **30" Wide Ductless Fume Hood** combines powerful air flow with a sturdy and compact enclosure for optimal respiratory and environmental protection. This fume hood helps protect the operator and environment from hazardous fumes or particulate generated from applications performed within the hood.

The **30" Ductless Fume Hood's** powerful fan directs fumes into the filter chamber containing high quality filter media selected depending on the application. Filter media choices include **HEPA** [up to 99.97% efficiency on particles 0.3 microns and larger], **ASHRAE** [up to 95% efficiency on particles 0.5 microns and larger, **Activated Carbon**, and **specialty-blended filtration** [i.e. acid gas, mercury, aldehyde, and ammonia].

Typical uses for the **30" Ductless Fume Hood** include chemical fumes, pharmaceutical compounding, soldering applications, dust removal, biological applications, and a variety of applications requiring fume and dust removal. The **30" Ductless Fume Hood** can provide an effective and economical solution for many commercial, medical, laboratory and industrial applications that require a small containment area. For additional visibility, a transparent version is available at extra cost.

Our line of ductless fume hoods is not intended for sterile compounding applications.





30" Ductless Fume Hood



30" DCH provides fume control for a various medical applications.

## **Product Features**

- Includes: variable speed controller (115V only), velometer, and LED light
- Optional: clear anti-static vinyl curtains, Magnehelic Gauge, hour counter, and spill tray
- **Standard Sizes:** 12", 18", 24", 30", 40", 50", 60", and 70" (custom sizes available)
- · Quiet operation
- · Long filter life
- · Simple, quick "no tool" filter change
- Portable and lightweight
- Reliable, low maintenance operation



## **Product Specifications**





Shown with outer dimensions. Dimensions are approximate.

CACE
JAJ
SENTRY AIR SYSTEMS, INC.

21221 FM 529 Rd Cypress, TX 77433 1.800.799.4609 sales@sentryair.com

www.sentryair.com

Hood Dimensions (D x W x H)  Outer: 24" x 32" x 25" Interior: 23.5" x 30" x 24" Access Area: 8" x 30"
(D x W x H) Interior: 23.5" x 30" x 24"
Access Area: 8" x 30"
<b>Sides:</b> 0.25" white acrylic
<b>Base:</b> 0.5" white high
density polyethylene
<b>Top:</b> 0.5" clear polycarbonate
Front Window: 0.25"
clear polycarbonate
Blower Dimensions 15.5"L x 13" W x 14.5" H
16 Ga. carbon steel with polyester epoxy coating
on all surfaces
Weight 100 lbs
Air Volume 350 CFM
(approximately) (varies with filter media)
Carbon Pre-filter + HEPA
Filter – 130 FPM high,
Inlet Velocities on 25 FPM low
High (approximately) Carbon Pre-filter + Heavy
Duty Carbon Filter –
115 FPM high, 20 FPM low
115/1/60, 2 amps
<b>Electrical</b> 220/1/50, 0.75 amps
8' grounded power cord
with NEMA 5-15P plug
Sound Level 61 to 66 dba (low to high)
Maximum Operating temperature
Temperature not to exceed 170° F
<b>HEPA</b> (up to 99.97%
efficiency on particles down to 0.3 microns)
,
<b>ASHRAE</b> (up to 95% efficiency on particles
Main Filter down to 0.5 microns)
Options Activated Carbon
Specialty-Blended Filter Media (i.e. acid gas,
mercury, aldehyde, and
ammonia)
Limited two-year warranty
from date of shinment on
Warranty defects due to materials
or workmanship.
<b>Patent #</b> 5,843,197