



HORIZONTAL LAMINAR FLOW HOOD WITH CONTROL CONSOLE

MODEL# SS-330-LFH



SENTRY AIR SYSTEMS **HORIZONTAL LAMINAR FLOW HOOD** offers protection to the operator from processes occurring within the hood. The horizontal laminar flow hood provides a negative pressure to the enclosure by filtering the contaminated air through HEPA and/or Carbon Filtration exhausting clean air back into the room.



- On/Off Switch
- Variable Speed Controller
- Flow Switch Indicator Light



- Air Velocity Alarm
- Hour Counter

CONFORMS TO CDCP & NIH STANDARDS FOR CLASS I BIOLOGICAL SAFETY CABINET*

*Per "Primary Containment for Biohazards" 3rd Edition Dept. of Health & Human Services & National Institutes of Health

HOOD DIMENSIONS	O.D.: 24"D x 32"W x 25"H I.D.: 23.5"D x 30"W x 24"H
HOOD MATERIALS	Sides: White Acrylic Base: White High Density Polyethylene Top: Clear Polycarbonate Front: Clear Polycarbonate
FAN/FILTER DIMS.	13" L x 13" W x 14" H
FAN/FILTER CHASSIS	16 ga. carbon steel w/ polyester epoxy coating on all surfaces
TOTAL WEIGHT	Approximately 100lbs.
MOTOR	Fractional Horsepower, TEFC
ELECTRICAL	115/1/60, 2 Amps
SOUND LEVEL	Approx. 55 dba @ 3' from hood inlet
APPROXIMATE INLET VELOCITIES	Pre-Filter and HEPA Filter- 130/25 FPM High/Low Pre-Filter & Carbon Filter- 115/20 FPM High/Low
MAX. OPERATING TEMP.	Not to exceed 170 degrees Fahrenheit
PATENTED DESIGN	U.S. PATENT #5,843,197

*Limited one year warranty from date of shipment on defects due to materials or workmanship.
All specifications are nominal and are subject to change without notice.
Not intended for explosive or highly combustible applications.

FEATURES:

Front Filter Access Door

- Enables easy access to filters

Audible Air Flow Alarm

- Indicates High/Low Air Velocities

Hour Counter

- Indicates filter and unit lifetime, displayed in hours

Variable Speed Controller

- No pre-set speeds, adjusts air flow

Magnetic Filter Maintenance Chart

- To record filter maintenance & facilitate reordering

Optional Minihelic Gauge

- Assists in monitoring HEPA filter condition

Ask about our... **Competitive Price Guarantee**

SENTRY AIR SYSTEMS, INC.



USA: 1.800.799.4609

Intl: 1.713.690.2153

Fax: 1.713.690.7872

www.sentryair.com