

## Ducted Fume Hood- 30" Wide



### Model # SS-330-E

### Product Specifications

#### HOOD DIMENSIONS

Outer Dimensions: 24" D x 32" W x 25" H  
Interior Dimensions: 23.5" D x 30" W x 24" H  
Access Area: 8" x 30"

#### HOOD MATERIAL

Sides: 0.250" White Acrylic  
Base: 0.500" White High Density Polyethylene  
Top: 0.500" Clear Polycarbonate  
Front: 0.250" Clear Polycarbonate

#### EXHAUST CONNECTION

6" O.D. Exhaust Collar  
Ducting & Fan Not Included with this Product

#### WEIGHT

Approx. 50 lbs.

#### MAX. OPERATING TEMPERATURE

Not to exceed 170 degrees Fahrenheit

#### WARRANTY

Limited two-year warranty from date of shipment on defects due to materials or workmanship.  
PATENT #5,843,197



Shown with outer dimensions. Measurements are approximate.



Disclaimer: The end user is ultimately responsible for conforming to all applicable exhaust ventilation regulations. Consult local, state and federal agencies regarding the suitability of your application prior to venting fumes to the environment. It is your ultimate responsibility to ensure the safety and suitability of your application. For more information, contact a local safety professional (certified industrial hygienist) for conformance and suitability.

### Product Features

- Fluorescent Light Included
- Custom Size & Modification Options
- Optional Spill Trays
- Custom Material Options
- Other Standard Sizes: 18", 24", 40", 50", 60", 70"
- Optional Anti-Static Vinyl Strip Curtains

**High-Quality Materials & Design.** Sentry Air's 30" Wide Ducted Exhaust Hood may be an excellent choice for applications that require a compact, exterior-ducting fume hood. This benchtop unit features unique hinged folding side panels, a high-density polyethylene base, a clear polycarbonate top, and a polycarbonate hinged sloping front lid. It is extremely compact and light weight. Typically, exhaust hoods may be used in chemical laboratories— especially those evaporating large amounts of volatiles organics, chemicals with low odor thresholds, heavy powder and capsule filling, and a variety of other processes that require exterior fume exhaust.

The 30" Wide Ducted Exhaust Hood can be an effective and economical tool for air quality enhancement in many commercial, chemical, medical, laboratory, and industrial applications.

**Customization Options.** This Ducted Fume Hood allows for several custom features:

- Material Options (e.g. metal, esd-safe, PVC)
- Size Options (width, height, depth, access area dimensions)
- Custom Cutouts (e.g. sink, machinery, microscopes)
- Dual-Entry Configuration
- Custom Illumination Options (e.g. Black Light)
- Accessories for Hanging Equipment

This unit includes a 6" Exhaust Collar and does not include ducting or fan. Ask Sentry Air Systems Applications Specialists for assistance in selecting a suitable fan. Sentry Air Systems also manufactures fume hoods that include a fan and filtration media (Ductless Fume Hood).



1-800-799-4609  
www.sentryair.com  
sales@sentryair.com

